

MAMMOGRAPHY GUIDELINE



SUBJECT: SPECIMEN RADIOGRAPH

1. Wire Localization Specimens will be done in the CONVENTIONAL MAG mode for Microcalcifications. Compressing the specimen is not required, unless the specimen is very large. The Technique is selected once you enter the Specimen Caption.
2. Wire Localization Specimens will be done in the CONVENTIONAL NON-MAG mode for Masses. Compressing the specimen is not required, unless the specimen is very large. The technique is selected once you enter the Specimen Caption.
3. Paraffin Block Specimens will be done in the CONVENTIONAL MAG mode for Microcalcifications. Compressing the specimen is not required, **even if a side view of the specimen is requested.** The technique is selected once you enter the Specimen Caption.