

DR. DARLINGTON

BARIUM ENEMA

SCOUT FILM: AP Abdomen 14x17

SUPPLIES: Sol-O-Pak Enema Bag filled to 2000cc
Enema Tip with cuff and balloon inflater

PROCEDURE: Set Fluoro technic , open study on TENS
Tip Patient after Dr. checks Scout
START patient supine
Operate Barium flow per Dr.

OVERHEADS: AP ABDOMEN
AP X-WISE FOR FLEXURES
LPO
LPO ANGLED OUTLET
LATERAL RECTUM
POST-EVAC