

CT Maxillofacial Guideline

Indications: Traumatic injury, Pain

Oral Contrast: None

IV Contrast: If indicated; Yes (See MRG Contrast Media Administration Policy)

Positioning	Supine, head first
FOV	Patient size
Axial Reconstruction: Thickness/ Interval/ Algorithm	(1) 2.5mm/2.5mm/Sharp (2) 2.5mm/2.5mm/Soft
Start Location	Below mandible
End Location	Through frontal sinus
Gantry Angle	None
Reformations: Thickness/ Interval/ Algorithm	(1) Coronal 2mm/2mm/Soft (2) Sagittal 2mm/2mm/Soft
Phase of Respiration	N/A