

DR. DARLINGTON

ACBE

SCOUT FILM: AP Abdomen 14x17

**SUPPLIES: empty Super XL Enema bag
1000cc Warmed Liquid Polibar Plus
Balloon inflator and bulb**

**PROCEDURE: Set Fluoro technic, open study on TIMS
Tip Patient after Dr. checks Scout
Start Patient Prone
Operate Barium flow per Dr.
Assist Dr. draining Patient upright**

**OVERHEADS: AP ABDOMEN
AP X-WISE for flexures
LPO
LPO ANGLED outlet
BOTH DECUBS
X-TABLE LATERAL (TIP OUT)**