

# DR HOFSTETTER

## ACBE

**SCOUT:** AP abdomen

**SUPPLIES:** Warmed POLIBAR (heat jug in a bucket in sink)

Large AIRCON Enema Bag, filled to 1500 cc's

INFLATOR and BULB attached to tip (check balloon)

**PROCEDURE:** Patient Tipped, START SUPINE

Tech controls Barium flow

table UPRIGHT, Tech opens valve, drains barium

Close valve, Pt SUPINE again, AIR introduced

**OVERHEADS:** **TECH gives more air between films**

AP- 14 x 17

LPO- 14 x 17

LPO ANGLED OUTLET- 14 x 17

BOTH DECUBS- 14 x 17

LATERAL RECTUM- pt. PRONE (crosstable - tip out)

POST EVAC- 14 x 17 if asked for