

DR. HOFSTETTER

BE

SCOUT: AP Abdomen 14 x 17

SUPPLIES: Sol -O-Pake Enema bag filled to 2000 cc's warm H₂O
Attach Flexi-cuff Enema tip, check balloon w/ inflator

PROCEDURE: Show SCOUT, if OK, TIP the Patient
START SUPINE, Tech controls Barium flow

OVERHEADS: AP- 14 X 17
LPO- 14 X 17
LPO ANGLED OUTLET- 14 x 17
LATERAL RECTUM- 10 X 12 OR 14 X 17
POST EVAC