

DR. DARLINGTON
SBFT

SUPPLIES: **MIX 2 cups ULTRA R BARIUM**

PROCEDURE: *Give Patient 1 cup ULTRA R BARIUM*
TAKE immediate film **AP STOMACH**
Give Patient 2nd ULTRA R cup to drink

FLUORO ROOM set up w/ PATIENT
When Ileocecal Valve is visualized
supine w/ paddle inflated for Rad.

OVERHEADS: **15 MIN**
30 MIN
CHECK Films w/ Dr. for timed intervals