

## **HOFSTETTER**

### **UGI**

#### **SUPPLIES:**

2 cups of ¼ cup POLIBAR (thick) Ba  
6 oz diluted POLIBAR (thin) Ba with straw  
CRYSTALS w/ 1Tbsp of water (2 sm cups)

#### **PROCEDURE:**

Scout patient and show Dr. Give Patient crystals in water and one cup of thick Barium. Lie Patient down and roll over 1 time. Take overheads. Start Patient UPRIGHT, Give patient second cup thick Ba. Dr. uses foot pedal. Patient SUPINE dr. takes spots. Patient then RAO, give thin Ba with straw.

#### **OVERHEAD IMAGES:**

AP  
RAO /LPO  
PA  
RIGHT LATERAL

#### **NO FAMILY IN ROOM**